

Professional
Powder Equipment
Manufacturer



Powder
Equipment



Milling
Technology



Powder
Materials

TENCAN

Product Brochure



ROLLER BALL MILL JAR

Stainless steel roll mill jar

Stainless steel roll mill jar widely used for grinding metal and non-metal materials. Customizable vacuum jar options available for different needs.



<https://www.labballmill.com/products/grinding-series/roller-ball-mill-jar/stainless-steel-roll-mill-jar.html>



Generated at 2026-07-28 15:35:01



Product Overview



Stainless steel roll mill jar widely used for grinding metal and non-metal materials.
Customizable vacuum jar options available for different needs.








Product Introduction

The roller stainless steel grinding jar is a widely used ball mill jar and is a commonly used grinding jar for grinding metal materials and non-metallic materials. Different models of vacuum tanks can also be customized according to customer requirements.






Technical parameters

model	Outer diameter(mm)	Inner diameter (mm)	Height including cover (mm)
500ml	F88	F82	144
1000ml	F106	Φ100	172
2000ml	F131	F124	206
3000ml	F158	F152	217
5000ml	F178	F172	269
10L	Φ243	Φ233	307
15L	F271	Φ263	347

Product Features



Accessories & Customization

Accessories

Grinding jars, heating elements, sample holders, control modules and other matching accessories can be selected according to the product configuration.

Customization

For voltage, capacity, chamber size, process temperature or application requirements, please contact TENCAN for a suitable configuration.