

Professional
Powder Equipment
Manufacturer



Powder
Equipment



Milling
Technology



Powder
Materials

TENCAN

Product Brochure



ROLLER BALL MILL JAR

cetamic roll mill jar

Ceramic mill jar: stable chemistry, wear resistant. Grinds chemicals, pharmaceuticals, rare earths, silicates.

<https://www.labballmill.com/products/grinding-series/roller-ball-mill-jar/cetamic-roll-mill-jar.html>

Generated at 2026-07-28 15:35:33







Product Overview

Ceramic mill jar: stable chemistry, wear resistant. Grinds chemicals, pharmaceuticals, rare earths, silicates.



Product Introduction

It has good chemical stability and wear resistance and is suitable for crushing various chemical products, pharmaceutical products, rare earth materials and various silicate materials.




Technical parameters

model	Diameter(mm)	Height including cover (mm)	Shoulder height (mm)
5L	Φ220	260	Around 220
10L	Φ260	320	Around 260
15L	Φ300	360	Around 280
20L	Φ300	400	Around 300

Product Features

1. The tank streamline design is perfectly unified with the ball mill.
2. The tank wall is super thick, wear-resistant and durable.
3. Economical and cost-effective.
4. Wide range of application.



Accessories & Customization

Accessories

Grinding jars, heating elements, sample holders, control modules and other matching accessories can be selected according to the product configuration.

Customization

For voltage, capacity, chamber size, process temperature or application requirements, please contact TENCAN for a suitable configuration.