

Professional
Powder Equipment
Manufacturer



Powder
Equipment



Milling
Technology



Powder
Materials

TENCAN

Product Brochure



ROLLER BALL MILL JAR

corundum roll mill jar

Corundum roll mill jar with extreme hardness, second only to diamond. Ideal for high-grade grinding, paired with zirconia balls for superior wear resistance.

<https://www.labballmill.com/products/grinding-series/roller-ball-mill-jar/corundum-roll-mill-jar.html>

Generated at 2026-07-28 15:37:02






Product Overview

Corundum roll mill jar with extreme hardness, second only to diamond. Ideal for high-grade grinding, paired with zirconia balls for superior wear resistance.



Product Introduction

Corundum, whose name originated from India, is a mineralogy name. The hardness of corundum is second only to diamond. Corundum is mainly used for high-grade abrasive materials. Corundum jars and zirconia balls are a classic wear-resistant combination. In addition, we need to remind our dear customers that although the corundum jar has strong hardness, it is also fragile. Please prevent it from falling to the ground and high-intensity impact during use to prevent it from breaking. Vacuum tanks of various materials and models can be customized according to customer requirements.



Technical parameters

model	Outer diameter(mm)	Excluding cover height (mm)	Height including cover (mm)	Reference weight (kg)
500ml	100	120	129	1.349
1L	142	150	160	2.757
1.5L	146	170	177	2.951
2L	152	189	197	4.039
3L	180	219	229	5.778
5L	203	257	267	7.968
10L	246	334	344	12.317

Product Features

High temperature resistance, corrosion resistance, high strength, etc.

Accessories & Customization

Accessories

Grinding jars, heating elements, sample holders, control modules and other matching accessories can be selected according to the product configuration.

Customization

For voltage, capacity, chamber size, process temperature or application requirements, please contact TENCAN for a suitable configuration.